

Planmeca Compact i Touch

New generation look & feel



PLANMECA






Step into the new era of usability

The new generation Planmeca Compact i Touch has a fully integrated touch panel for excellent usability and workflow. The dental unit use is made simple by displaying clear symbols and text in the large display. Both the dentist and assistant have good access to the panel and benefit from its excellent viewing angles.

Effortless instrument adjustments

The key advantage of the Planmeca Compact i Touch panel is a more straightforward use of instrument features. The dynamic user-interface displays all functionalities related to the selected instrument. All necessary adjustments are done effortlessly in one screen.



Dental unit general view



Micromotor selected

External infection control made easy

Smooth and uniform surfaces of Planmeca Compact i guarantee easy cleaning and impeccable external infection control, crucial for safe dental treatment. Powder-painted epoxy surfaces are extremely durable and cleaning agent resistant. All critical components can also be autoclaved and the bowl can be detached for easy cleaning.

Always clean water

The advanced water treatment solutions, with automated flushing cycles and Planmeca Waterline Cleaning System, guarantee that all water entering patient's mouth is always clean. It ensures safe dental treatment for the patient and safe working environment for the dental team.

Leading technology for clean suction tubes

Planmeca's innovative Suction Tube Cleaning System guarantees proper functioning of the dental unit suction system and safe treatment environment.

The system automatically dilutes and supplies the cleaning agent in the suction tubes. The daily cleaning procedure is fast, simple and free of extra accessories.



Extensive legroom

The side-lifting chair and small cuspidor base maximise the legroom for both the dentist and assistant. The floor is left free for stools and the dental team can work as close as possible to the patient. The lean and narrow backrest design also contributes to excellent access to the patient's mouth.

Practical upright sitting position

To guarantee fluent workflow, easy patient entrance and exit is supported with an automatic legrest. This feature also provides excellent patient comfort and is practical in doctor-patient consultation or prosthodontic treatments.

Compact in size and perfectly balanced

The small instrument console of Planmeca Compact i can be placed near patient for easy instrument reach. Perfectly balanced instrument arms guarantee ergonomic instruments use: only a tiny movement is necessary to activate the instrument. Instrument arms can also be adjusted depending on instrument weight.



Tomorrow's solutions today

Planmeca dental units are the first to provide innovative integration of digital dental unit and software. When connecting your dental unit to Planmeca Romexis software, all user activities (e.g. instruments use, cleaning routines, chair movements etc) are automatically recorded. The generated data can then be used for remote assistance, preventive maintenance and quality-control purposes. This solution, only available from Planmeca, guarantees a long lifespan for the dental unit.

Clear benefits with networked dental units

All Planmeca Compact i dental units can be networked and connected to Planmeca Romexis. A practical advantage is that dentists or hygienists can easily adjust and save their preferred automatic chair positions and instrument presets. The user-specific settings are instantly available from every unit of the clinic making it an ideal solution for large practices or dental universities with multiple users.

Planmeca Oy
 Asentajankatu 6 | 00880 Helsinki | Finland
 tel. +358 20 7795 500 | fax +358 20 7795 555
 sales@planmeca.com | www.planmeca.com